

TASKI Sani Des conc

Taski Sani Des Conc

Description

TASKI Sani Des conc is a cleaner and disinfectant for all washroom surfaces

Key properties

Versatile cleaner and disinfectant

Disinfectant, effective against a wide range of micro organisms

Fresh fragrance

TASKI Sani Des conc is part of the unique EXACT dosing and packaging concept for accurate dosage and cost control

Benefits

Removes common washroom soils

Reduces the risk of cross-contamination

Leaves a pleasant fragrance

Convenient and controlled dosing through the Exact system

Use instructions

Dosage:

Minimum dosage: 2 caps for 10L solution in bucket (0.2%). Increase concentration for heavy soiling.

Spray cleaning method: 1 cap in spray bottle with 500ml water (2%)

Spray Cleaning:

1 Spray the solution directly onto surface (or onto a cloth for smaller areas, and wipe).

2. On heavy soil leave for at least 5 minutes.

3. Clean surface with cloth or sponge

4. Rinse with clean water and allow to dry

Bucket Cleaning:

1. Apply the solution with a cloth/sponge or mop

2. Leave for at least 5 minutes

3. Rinse with clean water and allow to dry

Technical data

Appearance:

Clear, red liquid

Density at



20°C: 1.050

pH-value neat: 10.2 - 11.2

pH-value for both solution rates: 9.5 - 10.5

The above data is typical of normal production and should not be taken as a specification.

W9a



TASKI Sani Des conc

Safe handling and storage information

Safe handling:

Rinse and dry hands after use. The use of a skin reconditioning cream, such as **Soft Care Derasoft H9**, may be beneficial. For prolonged contact, protection for the skin may be necessary.

Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet.

For professional use only.

Store in original closed containers away from extremes of temperature.

Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EC Directives 73/404/EEC and 73/405/EEC and their subsequent amendments.

Microbiological data

TASKI Sani Des conc passed EN 1276 at a dilution of 1% in hard water (300ppm as CaCO₃) and in dirty conditions (0.3% bovine albumin), and contact time of 5 minutes.

TASKI Sani Des conc will pass an European suspension test with high soil level (1.0% bovine albumin) at a 1% concentration and contact time of 5 minutes.

Available pack sizes

Article code	Pack size	Pack type
7514525	1 litre	Dosing bottle

Diversey Ltd

Weston Favell Centre

Northampton NN3 8PD

Tel. 0800 525 525

www.diversey.com